


EPFL

Flying Cellulo for rehabilitation




Louis Jouret
Professor: P. Dillenbourg
Advisor: V. Borja

Goal of the Project



Study the feasibility of a haptic device to assist in arm movement for healthy aging and rehabilitation



Start the design of a first proof of concept



Market analysis

Hocoma is one of the leaders in movement rehabilitation.



Armeo®Power [1]

For early arm rehabilitation.
Completely assists the
patient's arm.

Armeo®Spring [2]

Provides arm weight
support.

Armeo®Senso [3]

Only sensor-based. For
moderate to mild
impairments

Solution: Flying drone


- solution that can be brought home
- can support the weight of a hand
- affordable

Requirements

- hand must be safe from propellers
- maximal load of 1.5kg
- no battery needed
- compatible with Playstation/Wii/Xbox communication protocol




Steps of the project




Fluid Simulations


1st iteration



4th iteration



final design



Control Algorithms

PID

- Easiest to design
- Difficult to tune
- Poor performance




no satisfactory result

EPFL


LQR

- Easiest to design
- Easy to tune
- Poor performance



MPC

- Difficult to design
- Easy to tune
- Best performance
- Very resource intensive




Test of the motors



Goal

- Check if motor fullfills requirements
- Find the relationship between the command and the thrust
- Find the maximal thrust for a given power

Results



Architecture




Test of the controller

Limitations:

- hazards of untested flying drone
- power supply not sufficient to have a flying drone


Solution:

- only study inclination → use of a balljoint



Test of the controller

Controller off



Controller on



Thanks

[1] “Armeo® power.”

[https://www.hocoma.com/solutions/armeo-power/.](https://www.hocoma.com/solutions/armeo-power/)

[2] “Armeo® spring.”

[https://www.hocoma.com/solutions/armeo-spring/.](https://www.hocoma.com/solutions/armeo-spring/)

[3] “Armeo® senso.”

[https://www.hocoma.com/solutions/armeo-senso/.](https://www.hocoma.com/solutions/armeo-senso/)

